

MONOAKCRYL

Synthetic Absorbable
POLYGLECAPRONE
SUTURES

OAK TREE[®]
SURGICAL



POLYGLECAPRONE SUTURES

MONOAKCRYL Polygecaprone monofilament sutures are rapid-absorbing synthetic monofilament sutures made from Poly(glycolide-co-caprolactone) and are available in both dyed and undyed versions. When dyed, only FDA-approved color additives, such as D&C Violet No. 2, are utilized. The undyed suture appears in a natural beige color.

MONOAKCRYL meets all the requirements specified in the European Pharmacopoeia (EP) for synthetic absorbable monofilament sutures. These sutures are exceptionally smooth, facilitating their passage through tissues when compared to braided sutures. They are infrequently utilized for percutaneous skin closure and are generally avoided in areas of high tension due to their characteristics.





















MONOAKCRYL exhibits low tissue reactivity and maintains approximately 50 to 60% of its in-vivo tensile strength at one week after placement. Complete absorption within tissues typically occurs over a period of around 90 to 110 days.

MONOAKCRYL

Polyglecaprone Monofilament

OAK TREE[®]

SURGICAL

Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
8mm						
 1/4	 50cm Violet	MA608FD		MA408FD		
	3/8 	50cm Violet	MB608D			
	3/8 	50cm Violet	MM608			
9mm						
 3/8	 50cm Violet	MB609				
	1/2 	30cm Violet	MN609	MN509		
	3/8 	50cm Violet		MB509		
11mm						
 3/8	 50cm Violet	ME611	ME511	ME411		
	 50cm Violet	ME611D				
	1/2 	50cm Violet	MF611	MF511		
	1/2 	30cm Violet	MN611	MN511	MN411	
	str. 	50cm Violet			MS411D	
13mm						
 3/8	 45cm Violet		ME513S			
	 50cm Violet	ME613	ME513	ME413		
	1/2 	50cm Violet	MF613	MF513		
	1/2 	70cm Violet	MN613	MN513		
	3/8 	60cm Violet			MM413	

MONOAKCRYL

Polyglecaprone Monofilament



Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
14mm 1/2 	30cm Violet	MN614	MN514			
	50cm Violet			MN414		
15mm 3/8 	50cm Violet	ME615	ME515	ME415	ME315	
	50cm Violet	ME615D		ME415D		
	1/2		MF515			
	3/8 	75cm Violet			MM415	
17mm 1/2 	50cm Violet		MN517S	MN417S	MN317S	MN217
	75cm Violet	MN617	MN517	MN417	MN317	
	5/8 	50cm Violet			MT417	
	3/8 	50cm Violet		ME517		
	75cm Violet			ME417	ME317	
	1/2	50cm Violet		MF517	MF417	
18mm 1/2 	50cm Violet	MN618	MN518	MN418	MN318	
	3/8 	50cm Violet			ME418	ME318
	1/2	50cm Violet		MF518	MF418	MF318
19mm 3/8 	50cm Violet		ME519	ME419	ME319	

MONOAKCRYL

Polyglecaprone Monofilament



Needle Type	Length & Color	Suture U.S.P (Metric)				
		5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)
19mm 3/8	50cm Violet	MM519	MM419	MM319		
	50cm Violet	MN519	MN419	MN319		
21mm 3/8	50cm Violet		ME421	ME321	ME221	
	50cm Violet	MF521	MF421	MF321	MF221	
	50cm Violet 75cm Violet	MN521	MN421	MN321	MN221	MN021
22mm 1/2	75cm Violet		MN422	MN322	MN222	
23mm 1/2	75cm Violet		MN423	MN323	MN223	
24mm 1/2	75cm Violet		MN424			
	50cm Violet			MN324		
25mm 3/8	75cm Violet		ME425	ME325	ME225	
	75cm Violet		MN425	MN325	MN225	MN025
	75cm Violet			MN325F		
	75cm Violet				MT225	

MONOAKCRYL











Polyglecaprone Monofilament



Needle Type	Length & Color	Suture U.S.P (Metric)				
		4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
26mm						
⊙ 5/8	75cm Violet		MN326F	MN226F		
⊙ 1/2	75cm Violet	MN426	MN326S	MN226S	MN026	
REMOVABLE	90cm Violet			MN226		
	45cm Violet		MN326C (8strands)			
▽ 3/8	75cm Violet		ME326	ME226	ME026	
⊕ 1/2	75cm Violet			MT226		
28mm						
⊙ 1/2	75cm Violet	MN428	MN328	MN228	MN028	
29mm						
▽ 3/8	75cm Violet		ME329			
1/2	75cm Violet		MF329	MF229		
⊙ 3/8	75cm Violet			MM229	MM029	
1/2	75cm Violet		MN329	MN229	MN029	MN0129









MONOAKCRYL
Polyglecaprone Monofilament



Needle Type	Length & Color	Suture U.S.P (Metric)			
		2-0 (3)	0 (3.5)	1 (4)	2 (5)
45mm  1/2 	75cm Violet	MN245	MN045	MN0145	MN0245
46mm  1/2 	75cm Violet	MT246		MT0146	
 1/2 	75cm Violet		MN046	MN0146	MN0246
 1/2 	75cm Violet		MF046	MF0146	
48mm  1/2 	90cm Violet	MN248		MN0148	

MONOAKCRYL
Polyglecaprone Monofilament




Needle Type	Length & Color	Suture U.S.P (Metric)			
		2-0 (3)	0 (3.5)	1 (4)	2 (5)
48mm  1/2 	75cm Violet				MF0248
49mm  1/2 	75cm Violet	MN249	MN049	MN0149	MN0249
50mm  str. 	75cm Violet	MS250			
51mm  1/2 	75cm Violet		MF051	MF0151	

MONOAKCRYL


Polyglecaprone Monofilament

OAK TREE
SURGICAL

Needle Type	Length & Color	Suture U.S.P (Metric)			
		2-0 (3)	0 (3,5)	1 (4)	2 (5)
<p>56mm</p> <p>⊙ 1/2</p> 	75cm Violet	MN256	MN056		
<p>63mm</p> <p>⊙ 3/8</p> 	75cm Violet		MU063	MU0163	
<p>66mm</p> <p>⊙ 3/8</p> 	75cm Violet	MM266	MM066		
<p>1/2</p> 	75cm Violet		MN066	MN0166	MN0266

MONOAKCRYL
Polyglecaprone Monofilament

OAK TREE[®]
SURGICAL













Needle Type	Length & Color	Suture U.S.P (Metric)			
		2-0 (3)	0 (3,5)	1 (4)	2 (5)
<div style="border: 1px solid black; padding: 2px; display: inline-block;">85mm</div> <input checked="" type="radio"/> 3/8 	75cm Violet			MU0185	

MONOAKCRYL

Polyglecaprone Monofilament

OAK TREE[®]











SURGICAL

Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
8mm  1/4 	45cm Undyed	MWA608F				
11mm  3/8 	45cm Undyed	MWE611	MWE511			
13mm  3/8 	45cm Undyed	MWE613	MWE513	MWE413		
16mm  3/8 	45cm Undyed		MWE516	MWE416		
17mm  1/2 	75cm Undyed				MWN317	
19mm  3/8 	45cm Undyed 75cm Undyed		MWE519	MWE419S MWE419	MWE319S MWE319	

MONOAKCRYL

Polyglactone Monofilament

OAK TREE
SURGICAL

Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
24mm  3/8 	45cm Undyed			MWE424		
	75cm Undyed			MWE424S	MWE324	MWE224
26mm  3/8 	75cm Undyed			MWE426	MWE326	
 1/2 	75cm Undyed				MWN326	
48mm  str. 	75cm Undyed			MWS448		
60mm  str. 	75cm Undyed			MWS460	MWS360	