



ISO 13485  
CERTIFIED COMPANY

# OAKSILK

BRAIDED SILK  
SUTURES

**OAK TREE**<sup>®</sup>  
SURGICAL



## BRAIDED SILK SUTURES

OAKSILK is a silk braided multifilament non absorbable suture.

**OAKSILK** is made by a unique developed process of degumming. This process removes extraneous SERICIN of the raw silk filament perfectly. This unique process allows for tighter, more compact braid which significantly improves suture quality. Silk filaments are braided on twisted core uniformly for greater tensile strength. It also ensure pliability for easy handling and knot security.

<b>Material</b>	Braided fibers from the cocoon of the silkworm
<b>Structure</b>	Braided
<b>Color</b>	Black, White
<b>Coating</b>	By silicone
<b>Range</b>	USP 8/0 to 5
<b>Absorption</b>	Non absorption
<b>Sterilization</b>	EO Gas Sterilization
<b>Advantage</b>	<ul style="list-style-type: none"> <li>Good and extended tissue support</li> <li>Excellent handling properties</li> <li>High tensile strength</li> <li>Excellent knot security</li> <li>Soft passage through the tissue with minimal sawing, tissue drag and traumatism</li> <li>Good knot tie down/adjustability</li> <li>Reduced capillary action</li> <li>Guaranteed tightness and product sterility</li> </ul>
<b>Application</b>	Ophthalmic surgery, Dental surgery, Ligation

**OAKSILK**  
Silk - Braided



Needle Type	Length & Color	Suture U.S.P (Metric)				
		8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)
6mm ▽ 3/8  FINE 220MIC	30cm Black	SK830A40D				
	50cm Black		SK750A41D			
	50cm White		SW750B17			
7mm ▽ 3/8	50cm Black		SK750C07			
8mm ▽ 3/8	50cm Black		SK750C09	SK650C38		
	50cm Black		SK750C09D	SK650C38D		
	50cm Black		SK750A44			
	50cm Black		SK750A44D			
9mm ▽ 3/8	50cm Black		SK750C10	SK650C41		
	50cm Black			SK650C41D		
11mm ▽ 3/8	50cm Black		SK750D04	SK650D29	SK550D57	
	50cm White			SW650D29		
	50cm Black			SK650D29D		
	50cm Black			SK650C49		
	50cm Black			SK650C49D		
13mm ▽ 3/8	50cm Black			SK650E14	SK550E42	SK450E63
	50cm Black					SK450D71
	75cm Black				SK575D62	
14mm ○ 1/2	30cm Black					SK430D72

**OAKSILK**  
Silk - Braided












**OAK TREE**  
SURGICAL

Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
<b>15mm</b>						
▽ 3/8	50cm Black	SK650F06	SK550F30	SK450F48	SK350F72	
○ 3/8	50cm Black		SK550F26	SK450F50		
	75cm Black			SK475F50		
▽ 3/8	50cm Black		SK550E47	SK450E68		
<b>17mm</b>						
○ 1/2	50cm Black		SK550I08	SK450I18	SK350I33	
	75cm Black				SK375F73	
	75cm Black				SK375H52	
▽ 3/8	50cm Black		SK550H10	SK450H27	SK350H53	
<b>18mm</b>						
▽ 3/8	50cm Black		SK550G22	SK450G44	SK350G69	
	50cm Black			SK450F57		
○ 1/2	50cm Black			SK450G43	SK350G67	
<b>19mm</b>						
○ 1/2	50cm Black			SK450E73	SK350E88	
<b>21mm</b>						
▽ 3/8	50cm Black			SK450I26	SK350I36	
▽ 3/8	50cm Black			SK450H32	SK350H61	SK250H77
○ 3/8	75cm Black				SK375H58	
▽ 1/2	50cm Black				SK350J22	
○ 1/2	75cm Black			SK475J15	SK375J23	

# OAKSILK

Silk - Braided

**OAK TREE**  
SURGICAL

Needle Type	Length & Color	Suture U.S.P (Metric)				
		4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
<b>23mm</b> ○ 1/2 	75cm Black		SK375H63			
<b>24mm</b> ○ 1/2 	50cm Black		SK350G77			
<b>25mm</b> ▽ 3/8 	50cm Black		SK350J25	SK250J36		
○ 1/2 	75cm Black		SK375J24	SK275J34		
○ 1/2 	75cm Black		SK375I41	SK275I52		
○ 1/2  8 STREND/PACK	75cm Black 45cm Black		SK375G56	SK275H84		
<b>28mm</b> ○ 1/2 	75cm Black		SK375I43			
<b>29mm</b> ▽ 3/8 	75cm Black		SK375K11	SK275K16	SK175K26	
○ 1/2 	75cm Black		SK375K09			
<b>33mm</b> ▽ 3/8 	75cm Black			SK275K20	SK175K28	
○ 1/2 	75cm Black			SK275L13	SK175L23	SK0175L32

# OAKSILK

Silk - Braided



Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
41mm ▽ 3/8 	75cm Black				SK0175N47	

## LIGATURE AND NON-NEEDLE SUTURES

### ■ BRAIDED BLACK SILK

U.S.P (Metric)	BLACK SILK						
	180cm × 1	150cm × 1	110cm × 7	75cm × 10	60cm × 13	45cm × 17	45cm × 10
6-0 (0.7)				SK - 67510		SK - 64517	SK - 64510
5-0 (1)			SK - 51107	SK - 57510		SK - 54517	SK - 64510
4-0 (1.5)		SK - 41501	SK - 41107	SK - 47510	SK - 46013	SK - 44517	SK - 44510
3-0 (2)	SK - 31801	SK - 31501	SK - 31107	SK - 37510	SK - 36013	SK - 34517	SK - 34510
2-0 (3)	SK - 21801	SK - 21501	SK - 21107	SK - 27510	SK - 26013	SK - 24517	SK - 24510
0 (3.5)		SK - 11501	SK - 11107	SK - 17510	SK - 16013	SK - 14517	SK - 14510
1 (4)	SK - 011801	SK - 011501		SK - 017510	SK - 016013	SK - 014517	SK 014510
2 (5)				SK - 027510		SK - 024517	SK - 024510

### ■ BRAIDED BLACK SILK

U.S.P (Metric)	WHITE SILK				
	75cm × 10	60cm × 13	45cm × 10	40cm × 10	40cm × 5
4-0 (1.5)		SW - 46013	SW - 44510	SW - 44010	SW - 4405
3-0 (2)	SW - 37510	SW - 36013	SW - 34510	SW - 34010	SW - 3405
2-0 (3)	SW - 27510	SW - 26013	SW - 24510	SW - 24010	SW - 2405
0 (3.5)	SW - 17510	SW - 16013	SW - 14510	SW - 14010	SW - 1405
1 (4)	SW - 017510		SW - 014510	SW - 014010	SW - 01405
2 (5)				SW - 024010	SW - 02405





**OAK TREE<sup>®</sup>**  
SURGICAL